

## Nano-Clear - Accelerated Weathering Test - Using Super Xenon WOM

Sample Details	ΔE Values (Super Xenon WOM)							20° Gloss retention (%)		
	100 hours	200 hours	300 hours	400 hours	500 hours	600 hours	700 hours	500 hours	600 hours	700 hours
Nano-Clear	0.99 to 1.00	1.00 to 1.06	0.95 to 0.96	0.81 to 0.83	0.67 to 0.71	0.57 to 0.59	0.48 to 0.53	100%	102%	99%

